
öko 2000

batch no.:-----

programme no. 1

1: cleaning

beginning: 01.01.70 03:27

duration: 20 min.

temperature: 45 °C

DOS1: neutral

Komponente A 88 ml

2: cleaning

beginning: 01.01.70 03:57

duration: 20 min.

temperature: 40 °C

with detergent powder

3: rinsing

beginning: 01.01.70 04:25

duration: 4 min.

temperature: 40 °C

ec-value: none

4: rinsing

beginning: 01.01.70 04:35

duration: 4 min.

temperature: 30 °C

ec-value: 193 µS

5: rinsing

beginning: 01.01.70 04:43

duration: 4 min.

temperature: 30 °C

ec-value: 37 µS

6: rinsing

beginning: 01.01.70 04:51

duration: 4 min.

temperature: 30 °C

ec-value: 10 µS

7: drying

beginning: 01.01.70 05:00

run after: 20 min.

temperature: 70 °C

dry grade: 5

ending: 01.01.70 05:39

failure: none

name:-----