

C67x FastRTS vs. RTS Benchmarking

Single Precision				
Function	RTS6740**	FastRTS	FastRTS (Intrinsics)	FastRTS (Inlining) Pipelining w/128 Calls
atan	651	91	101	19
atan2	1267	89	100	26
cos	167	74	94	18
div	540	37	38	3
exp	184	82	95	18
exp2	190	82	143	22
exp10	190	82	143	22
log	128	74	75	13
log2	153	75	79	17
log10	153	75	79	17
pow	937	181	177	82*
recip	542	32	42	3
rsqrt	165	44	78	18
sin	142	69	82	17
sqrt	161	49	70	17

Double Precision				
Function	RTS6740**	FastRTS	FastRTS (Intrinsics)	FastRTS (Inlining) Pipelining w/128 Calls
atan	1915	320	356	322
atan2	2016	453	479	421
cos	327	155	190	103
div	853	99	118	61
exp	1136	229	287	96
exp2	1083	243	297	96
exp10	1157	243	297	96
log	1871	302	359	326
log2	1887	305	370	337
log10	1891	305	370	337
pow	1570	552	587	583
recip	856	90	108	61
rsqrt	325	120	139	45
sin	309	150	188	104
sqrt	337	130	151	81*

Notes: The benchmarking was performed using the RTS, fastRTS, fastRTS w/intrinsics and fastRTS w/inlining of the different math functions. The inlining was performed to allow for pipelining in a loop iterated 128 times. The average clock cycles was then determined for the respective fastRTS (Inlining) column. The functions atan2 and pow require 2 arguments. *Special case modified to optimize. ** Average clocks measured may vary depending on the range of the data

